

# CDRWineLab®

Analysis system  
for vinegar  
quality control



## CDR WineLab® system

CDR WineLab® consists of a thermostatically controlled analyser with photometric technology using LED emitters and kits of reagents that are pre-filled into vials and ready to use.



1

### Take the sample

to be analysed using the pipettes supplied with the system.



2

### Place the sample

in the test tube containing the pre-filled reagent.



3

### Insert the test tube

into the reading cell to obtain the analysis result.

CDR WineLab®, the wine analysis system, can also be configured to carry out analyses on **wine vinegar**. The system thus becomes a valuable aid for **quality control** in **vinegar production**.



## Reduced analysis times

With CDR WineLab® you are finally free to carry out the analyses independently, in your own wine cellar, quickly and easily, without having to rely on an external laboratory. In fact, it is possible to analyse **16 samples simultaneously** and constantly monitor the production process, obtaining specific and precise answers in a few minutes.



## Easy to use

The system has been designed so that it can be used **not only in the laboratory, but also on the production line for real-time results**, by personnel without specific technical training.

**The analysis methods, shown on the display, are simpler than traditional methods and can be performed in just a few steps.**

If required, the HELP function will guide the operator step by step through the procedure. The result is automatically calculated, displayed and printed out.




## Reliable

CDR WineLab® **guarantees high sensitivity, a wide measuring range and excellent repeatability** of the results thanks to the innovative photometric technology using LED light sources and fixed wavelengths ranging from the ultraviolet to the visible spectrum (with a range of 0 to 6 optical density). **The analysis results are correlated with those of the reference methods.**

Pre-filled and disposable reagents are packaged in bags of 10 tests, developed and produced by the CDR research laboratories.





## Analyses of vinegar with CDR WineLab®

	TEST	Measuring range	Resolution	Repeatability	Test time	
Vinegar	Total acidity	2.50-12.00 g/100 ml of acetic acid	0.01 g/100 ml of acetic acid	0.30 g/100 ml of acetic acid	1 min	
	Alcohol content	0.10-3.50 % vol.	0.01 % vol	0.03 % vol	11 mins	
		3.0-7.0 % vol.	0.1 % vol	0.1 % vol	11 mins	
	Total SO <sub>2</sub> on white vinegar	15-250 mg/l	1 mg/l	4 mg/l	1 min	
	Total SO <sub>2</sub> on red vinegar	15-250 mg/l	1 mg/l	6 mg/l	1 min	
	Sugars	30-600 g/l of glucose	1 g/l of glucose	4 g/l of glucose	6 mins	
		400-1000 g/l of glucose	1 g/l of glucose	7 g/l of glucose		
	Chlorides	50-400 mg/dl NaCl	1 mg/dl NaCl	15 mg/dl NaCl	5 mins	
	Colour (intensity and tonality)	0.000-40.000 Intensity	0.001	0.100	1 min	
			0.000 - ∞ Tonality	0.001		0.050
			1.0 - 5.000 - 420 nm reading	0.001		0.072
			1.0 - 5.000 - 520 nm reading	0.001		0.072
1.0 - 5.000 - 620 nm reading			0.001	0.072		
15. - 100.00 Intensity		0.01	0.10	1 min		
		1.0 - ∞ Tonality	0.001		0.050	
		1.0 -5.000 - 420 nm reading	0.01		0.25	
		1.0 - 5.000 - 520 nm reading	0.01		0.25	
		1.0 - 5.000 - 620 nm reading	0.01		0.25	

Faster and more intuitive than traditional instruments, **CDR WineLab®** can be used to improve and accelerate standard quality controls in **vinegar production**.

# CDRWineLab®



# CDRWineLab® Jr



<b>Analyses</b>	Complete analysis panel	Customisable configuration
<b>Samples that can be analysed simultaneously</b>	16	3
<b>Multitasking Mode</b>	Yes	No
<b>Calibration</b>	Pre-calibrated No periodic calibration is necessary	Pre-calibrated No periodic calibration is necessary
<b>Maintenance costs</b>	No	No
<b>Storage of results</b>	Sufficient internal memory for storing thousands of analysis results in CVS and XML files compatible with all database formats (e.g., XLS, SQL)	Sufficient internal memory for storing thousands of analysis results in CVS and XML files compatible with all database formats (e.g., XLS, SQL)
<b>Photometric module</b>	Up to 6 wavelengths in 4 reading cells	Up to 6 wavelengths in 4 reading cells
<b>Incubation module</b>	37 ° C thermostated block with 16 positions	37°C thermostated reading block with 3 positions with incubation function
<b>Connection with barcode and QR code scanners</b>	Yes, via Bluetooth	No
<b>Display</b>	5.7" TFT colour LCD with touch screen	4.3" TFT colour LCD with touch screen
<b>Connectivity</b>	1 USB port type B for transferring the performed analysis database, configuration and software update, PC connection 1 USB port type A for technical service and computer connection 1 Ethernet port (LAN) for connection to intranet Bluetooth 4.0	1 USB port type B for transferring the performed analysis database, configuration and software update, PC connection  Bluetooth 2.1
<b>Printer</b>	80 mm wide printer with integrated graphics	Wireless connection for external printer
<b>Dimensions and weight</b>	32 x 29.5 x 13 cm (W x D x H) 2.80 kg	15 x 22 x 8,3 cm (W x D x H) 0,80 Kg
<b>Power supply</b>	24 V	24 V or optional lithium-ion battery

rev.7.0

## CDRWineLab®

CDR WineLab®, system of  FOODLAB® line, is a trademark of CDR S.r.l. Phone. +39.055.871431 • Fax +39.055.8714322  
Via degli Artigiani, 6 • 50055 Ginestra Fiorentina • Florence - ITALY • info@cdrwinelab.com - www.cdrwinelab.com

